Venus® Bulk Fill

Produktbeschreibung

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Venus® Bulk Fill

Product description

Venus Bulk Fill is a posterior flowable, light-curing, radiopaque nano-hybrid composite. It is used as a base in Class I and II restorations and can be placed in 4 mm increments. Venus Bulk Fill's flow and demand characteristics enable the material to adapt to the cavity walls. The universal shade allows curing of 4 mm layer thickness (base fill) with low shrinkage stress.

A universal restorative material is used in conjunction with Venus Bulk Fill as a facial/oral/occlusal enamel replacement. Venus Bulk Fill is compatible with methacrylate-based restorative materials designed for Class I and II occlusal restorations (e.g. Venus Diamond® see Application).

Composition

Venus Bulk Fill is based on multifunctional methacrylate monomers (UDMA, EBADMA) and contains approximately 65 % w/w and 38 % vol inorganic fillers, such as Ba-Al-F silicate glass, YbF₃, and SiO₂. The filler particle size is between 0.02 µm and 5 µm.

Color

Universal shade

Indications

- Base in Class I and II direct restorations (up to 4 mm)
- Cavity lining under direct restorative materials in Class I and II

Contraindications

Do not use in cases of a known allergy to any of the constituents of Venus Bulk Fill.

Application to the dentin less than 1 mm to the pulp

Side effects

Allergies to the product or its constituents may occur in isolated cases. If an allergy is suspected, request a list of ingredients from the manufacturer.

Warnings / Precautions

Avoid contact with skin or eyes. May cause sensitisation in case of skin contact. May irritate eyes. In case of contact with skin, wash skin with soap and water. In case of contact with eyes, rinse with plenty of water and seek medical attention.

PLTs are intended for single use. To avoid cross-contamination between patients, PLTs must be discarded after a single use!

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For deep cavities with only a minimal layer of dentin over the pulp or those which have already had a small area of pulp exposed, the pulp must be indirectly or directly capped locally using calcium hydroxide. The calcium hydroxide should be covered with a glass ionomer cement before Venus Bulk Fill is applied to the cavity.

Prepared enamel and dentin surfaces on which Venus Bulk Fill will be adhesively luted must not come in contact with preparations containing eug-enol, as these could inhibit optimal polymerization of the composite.

Deviations from the instructions for use are the sole responsibility of the dentist. The product must be used solely in accordance with the instructions for use.

Application

NOTE: The shade selection of the final composite should be selected using the shade guide before inserting the dental dam as teeth noticeably change color when dried.

- Remove caries, residual amalgam or restorative material from the cavity. Then rinse, clean and dry the cavity but do not desiccate the tooth structure.
- Use a dental dam (e.g. Ivory® Rubber Dam) to isolate the treatment area from contamination. Refer to manufacturer's directions for placement.
- Shape the proximal contacts with the aid of a matrix system. The use of a matrix band with a wedge serves two purposes: compressing matrix band against the tooth structure creates a tighter seal, preventing gingival overhang and it provides slight separation between the teeth, compensating for the thickness of the matrix band, therefore helping create tight contact between the finished restoration and the adjacent tooth. Therefore, a thin matrix band is recommended. The matrix band should be burnished, or shaped, to achieve a contour facilitating its proper adaptation to the tooth structure. The matrix band should also be shaped with a suitable instrument while applying the composite, because Venus Bulk Fill does not exert pressure during application.
- An appropriate adhesive system (e.g., iBond® Self Etch, iBond® Total Etch) should be applied to the enamel and dentin according to the manufacturer's instructions before applying the filling material. It is the responsibility of the dental practitioner if another dentin/enamel adhesive product is utilized. Adhesive-treated cavities must not be contaminated with moisture, blood or saliva.
- Begin dispensing Venus Bulk Fill directly into the cavity preparation at the deepest point, keeping tip submerged. When applying, ensure that it flows slowly and evenly to prevent the formation of air bubbles. Immerse the PLT tip into the material and fill the cavity from bottom to top. If necessary, use a probe or a micro-brush to remove any air bubbles and excess material before polymerization. The distance to the occlusal margin of the cavity must be at least 2 mm.
- Alternative application would be to use Venus Bulk Fill in a thin layer as a liner on exposed dentin.

Venus Bulk Fill is designed to be cured in increments up to 4mm in depth. Polymerize Venus Bulk Fill for 20 seconds, on each of the restoration surface with a curing light (e.g., Translux® Power Blue) in accordance with the table below (light output: > 550 mW/cm²).

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Stand: 05.2010

When using plasma arc curing lights, the polymerization time can usually be reduced by 25 %. When using LED (light emitting diode) curing lights, follow the polymerization times specified for curing with halogen units, provided direct contact between the light source and the composite can be ensured. If the light intensity falls below 550 mW/cm², increasing the time by 50 % is recommended. It must be noted that polymerization lamps can vary significantly in intensity and emitted spectrum. All polymerization times specified are therefore minimum times. In case of doubt, contact the lamp manufacturer for more information and check the functionality of the curing light in vitro before using it on patients.

During polymerization, a dispersion layer forms on the surface. This layer must not be touched or removed as it ensures optimum application of further composite layers.

Complete the restoration with a universal or posterior tooth composite (e.g., Venus Diamond) without delay. In general, the choice of a Universal Shade is the best for an aesthetic result.

If it is necessary to contour or remove cured Venus Bulk Fill, use rotary instruments (e.g., carbide instruments) for trimming as directed by the manufacturer.

Instructions

During long restorative procedures, the surgery light should be temporarily directed away from the treatment area to prevent premature initiation of the polymerization reaction.

The polymerization time specified above is a minimum time and should never be reduced to allow an optimal polymerization rate and the associated physical material properties to be achieved.

Please inspect for damage before using for the first time. Do not use damaged products.

Venus Bulk Fill should be used at room temperature (23 °C/73 °F).

Further precautions

Only to be used by dentists and for its intended use.

Storage

Store below 25 °/77 °. Keep out of direct sunlight.

Do not use the material after the expiry date.

Expiry date ☞ and batch number ☞ See note on syringe/packageing.

Please quote batch number and expiry date in all correspondence about the product.

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